

## Purifying Amplified Products

The C-WTA products should be purified to remove residual primers and nucleotides which may interfere with downstream processes, such as labeling reactions.

The C-WTA product should be purified on one QIAquick column as instructed in the QIAquick® PCR Purification Kit (QIAGEN, Cat # 28104 or 28106), with the following modifications to the standard protocol:

1. Add 375 µL of Buffer PB and 10 µL of 3M sodium acetate (pH 5.2) to the 75 µL product and mix by vortexing before loading the sample on a column.
2. Elute purified product in 50 µL of TE Buffer (pH 8.0), allowing column to stand for 2 minutes at room temperature after adding TE before centrifuging.

The purified, amplified product may be stored at -20°C.

## Quantifying Amplified Products

UV absorbance (A260) should be used to quantify purified products, using the conversion of 1 OD = 50 µg/mL. Approximately 10-15 µg of amplified product should be generated from 10 ng of Universal Human Reference RNA (Agilent, Cat # 740000) sample.

SPE-LBL-0250.00



**Product Number: RC20050**  
**Amount: 50 reactions**  
**Storage: -20°C**  
**FOR RESEARCH USE ONLY**

## TransPLEX® C-WTA Kit

**For Whole Transcriptome Amplification of Clinical RNA Samples**

### Kit Components

Component Name (Part Number)	Cap Color	Volume
C-WTA Synthesis Buffer (RC20050-01)	Green	325 µL
C-WTA Stabilization Solution (RC20050-02)	Yellow	140 µL
C-WTA Synthesis Enzyme (RC20050-03)	Blue	80 µL
C-WTA Amplification Buffer (RC20050-04)	(White)	825 µL
C-WTA Amplification Enzyme (RC20050-05)	Red	65 µL
Nuclease-Free Water (RC20050-06)	(Clear)	3 x 1.8 mL
User Manual (RC20050-07, v1)		

### Storage and Handling

Store the TransPLEX C-WTA Kit at -20°C. Transfer C-WTA Synthesis Enzyme and C-WTA Amplification Enzyme to ice and briefly centrifuge just before use. Thaw other components on ice and briefly vortex and centrifuge prior to use.

### User-Supplied Materials

- Thermal cycler (Real-time instrument recommended)
- PCR tubes or 96-well PCR plate
- PCR plate seals
- Low-binding barrier tips
- QIAquick® PCR Purification Kit (QIAGEN, Cat # 28104 or 28106 )

TransPLEX® is covered by US Patent 7,655,791 and related US and foreign patents.

TransPLEX® is a registered trademark of Rubicon Genomics, Inc.  
QIAquick® is a registered trademark of QIAGEN.


Rubicon Genomics, Inc., 4355 Varsity Drive, Ann Arbor, MI 48108 USA  
Tel: 734-677-4845 ▪ Internet: [www.rubicongenomics.com](http://www.rubicongenomics.com)

## Introduction

The TransPLEX C-WTA Kit allows amplification of total RNA to over 5 micrograms of amplified DNA products suitable for both microarray and qPCR analyses in 5 hours.

TransPLEX C-WTA has two steps, library synthesis and amplification. To synthesize the library, sample RNA is incubated with a reverse transcriptase and non-self-complementary primers comprised of a quasi-random 3' end and a universal 5' end. When annealed primers are extended by polymerase, displaced single strands are generated which become new templates for primer annealing and extension. This process creates a TransPLEX library comprised of random, overlapping fragments flanked by a universal end sequence. Universal-primer PCR is then used to amplify the TransPLEX library and produce C-WTA products.

## TransPLEX C-WTA Kit Sample Processing Flow



## TransPLEX C-WTA Kit Protocol (Library Synthesis)

1. Add Nuclease-Free Water to total RNA (10 ng of high quality or 25-300 ng of FFPE) to achieve a total sample volume of 16.5 µL in a PCR tube or well.
2. Combine the following Library Synthesis Cocktail components and mix well.

Library Synthesis Cocktail Component	Volume Per Sample
C-WTA Synthesis Buffer	5 µL
C-WTA Stabilization Solution	2.5 µL
<b>Total Volume</b>	<b>7.5 µL</b>

3. Add 7.5 µL of freshly prepared Library Synthesis Cocktail to the RNA sample and mix by pipet.
4. Incubate sample in a thermal cycler as follows:

1 cycle	70°C	5 min
1 cycle	4°C	Hold

5. Briefly centrifuge sample to collect liquid at bottom of tube/well and place sample on ice.
6. Add 1 µL of C-WTA Synthesis Enzyme to sample for a total of volume of 25 µL and mix by pipet.
7. Incubate sample in a thermal cycler as follows to produce a TransPLEX Library:

1 cycle	24°C	15 min
1 cycle	42°C	2 hour
1 cycle	95°C	5 min
1 cycle	4°C	Hold

8. Briefly centrifuge the Library and transfer a single 5 µL aliquot of the Library to a new tube/well to prepare for the following Amplification procedure.

## TransPLEX C-WTA Kit Protocol (Library Amplification)

9. Combine the following Amplification Cocktail components and mix well.

Amplification Cocktail Component	Volume Per Library Aliquot
Nuclease-Free Water	54 µL
C-WTA Amplification Buffer	15 µL
C-WTA Amplification Enzyme	1 µL
<b>Total Volume</b>	<b>70 µL</b>

10. Add 70 µL of freshly prepared Amplification Cocktail to the 5 µL Library aliquot (prepared in Step 8) and mix by pipet.
11. Amplify sample according to thermal cycler program below:

1 cycle	95°C	2 min
22 cycles	95°C	20 sec
	65°C	5 min
1 cycle	4°C	Hold

12. Briefly centrifuge and immediately store the amplified product at -20°C or proceed to purifying the amplified products.