

Protocol: Target Enrichment

Exome Capture of ThruPLEX[®] Libraries with Agilent SureSelect[®]QXT Target Enrichment System

For use with all ThruPLEX DNA-seq Kits and ThruPLEX Plasma-seq Kits

Required Reagents

ThruPLEX[®] library preparation kit (one of following required):

- ThruPLEX DNA-seq Kit (Rubicon Genomics)
 - 12 rxns, 12 single indexes (CAT. NO. R400429)
 - 48 rxns, 12 single indexes (CAT. NO. R400428)
 - 48 rxns, 48 single indexes (CAT. NO. R400427)
 - 48 rxns, 48 dual indexes (CAT. NO. R400406)
 - 96 rxns, 96 dual indexes (CAT. NO. R400407)
- ThruPLEX Plasma-seq Kit (Rubicon Genomics)
 - 12 rxns, 12 single indexes (CAT. NO. R400490)
 - 48 rxns, 48 single indexes (CAT. NO. R400491)
 - 96 rxns, 96 dual indexes (CAT. NO. R400492)

Blocking oligos (both required):

- xGen[®] Universal Blocking Oligo – TS HT-i5 (Integrated DNA Technologies)
- xGen Universal Blocking Oligo – TS HT-i7 (Integrated DNA Technologies)

Primers (both required):

- Illumina[®] P5 Primer: AATGATACGGCGACCACCGA
- Illumina P7 Primer: CAAGCAGAAGACGGCATACGA

SureSelect[®]QXT reagents:

- Refer to Required Reagents on page 13 in the Agilent SureSelectQXT Protocol (Version B.2, October 2014)

Required Equipment

Refer to Required Equipment on page 15 in the Agilent SureSelectQXT Protocol (Version B.2, October 2014)

Important Notes

When integrating ThruPLEX with the Agilent SureSelectQXT2 library capture system, all components of the SureSelectQXT Reagent Kit are used except the following:

- SureSelectQXT Buffer
- SureSelectQXT Enzyme Mix ILM
- DMSO
- SureSelectQXT Read Primer 1

- SureSelectQXT Read Primer 2
- SureSelectQXT Index Read Primer
- SureSelectQXT P7 dual indexing primers
- SureSelectQXT P5 dual indexing primers
- SureSelectQXT Stop Solution
- SureSelectQXT Primer Mix

Contact Agilent to order a SureSelectQXT Reagent Kit without the SureSelectQXT Library Prep Kit, ILM, Box 2.

CAUTION: This custom kit may require the following item for the post-capture amplification step:

- Herculase II Fusion DNA Polymerase with dNTPs (Agilent Technologies, CAT. NO. 600677 or 600679)

ThruPLEX Library Preparation

1. Prepare ThruPLEX libraries
 - Follow Section D in the ThruPLEX DNA-seq Kit or ThruPLEX Plasma-seq Kit Instruction Manual.
2. Perform Library Purification by AMPure XP beads
 - Follow Section E.V. in the ThruPLEX DNA-seq Kit or ThruPLEX Plasma-seq Kit Instruction Manual.
 - **CAUTION:** For the final elution, DNA **must** be eluted by resuspending the beads in 30 μ L of PCR grade water, **not** TE buffer.

Capture of ThruPLEX Libraries

1. Resuspend xGen Universal Blocking Oligos to 1 μ L per reaction (or 1 nmole/ μ L) in nuclease-free water.
2. Using a narrow gauge needle, poke hole(s) in the lid of each tube containing a ThruPLEX library to be used for capture.
3. Concentrate the ThruPLEX library using a vacuum concentrator held at $\leq 45^{\circ}\text{C}$ to reduce the volume in the tube to $< 10 \mu\text{L}$. Do not completely dry the mixture.
4. Bring the volume to 10 μ L with nuclease-free water.

5. To each resuspended library, add:
 - 1 μ L of xGen Universal Blocking Oligo – TS HT-i5
 - 1 μ L of xGen Universal Blocking Oligo – TS HT-i7
6. Follow procedures in the Agilent SureSelectQXT Protocol (Version B.2, October 2014) starting at Chapter 3, Step 2 through the end of Chapter 4, Step 5 with the following modification:

Chapter 4, Step 1. Amplify the Captured Libraries

- Modify the Post-Capture PCR Reaction Mix (Table 18, page 47) to the following:

Reagent	Volume for 1 rxn
Nuclease-free water*	10.5 μ L
5x Herculase Rxn Buffer	10.0 μ L
100 mM dNTP Mix	0.5 μ L
Herculase II Fusion DNA Polymerase	1.0 μ L
10 μ M Illumina P5 Primer	2.5 μ L
10 μ M Illumina P7 Primer	2.5 μ L
Total	27.0 μL

- ThruPLEX libraries are already indexed, so do **not** use the SureSelectQXT indexing primers.

Note: This protocol was developed using the SureSelectXT Human All Exon v5 Capture Library.

ThruPLEX® is a registered trademark of Rubicon Genomics, Inc., SureSelect® is a registered trademark of Agilent Technologies, Inc., xGen® is a registered trademark of Integrated DNA Technologies, Inc., Agencourt® AMPure® are registered trademarks of Beckman Coulter, Inc.

ThruPLEX® Kits are intended for Research Use Only. It may not be used for any other purposes including, but not limited to, use in diagnostics, forensics, therapeutics, or in humans. ThruPLEX may not be transferred to third parties, resold, modified for resale or used to manufacture commercial products without prior written approval of Rubicon Genomics, Inc.

ThruPLEX DNA-seq and ThruPLEX Plasma-seq are protected by U.S. Patents 7,803,550; 8,071,312; 8,399,199; 8,728,737 and corresponding foreign patents. Additional patents pending.

